





In-situ high-pressure electrochemistry



Well-defined force application

- Uniaxial force application
- Load-bearing parts made from hard-metal
- Pressures up to 1.7 GPa (Ø 6 mm, 50 kN)
- Force application via high-precision servo-motor instead of hydraulics




Ideal for electrochemical *in-situ* characterization

- Plate-plate geometry
- > Connectable with all potentiostats
- > Sandwich-type material stacks possible
- > Electrical contact via plunger and bottom plate
- > Chemically inert & electrically insulated inner sleeve



Precisely controlled pressure and temperature conditions





Active force regulation

- Applied force is constantly measured and adjusted as long as required
- > Constant within ±0.1 kN
- Force setpoints can be applied manually, via procedure or software interface

Accurate temperature control

- Possible setups: passive/ electrical heating/ fluid cooling & heating*
- Enables temperature-dependent analyses from -40 °C* to +250 °C** (±0.1 °C)
- High quality Pt-100 temperature sensor







CompreDrive H CompreDrive HC* *external circulator required **depending on sleeve material

rhd instruments


The CompreCell – Handling of sensitive samples


Easy to use, air-tight setup

rhd instruments

flexible cell solutions

- > Handling of air-/moisture-sensitive samples
- > Easy, glovebox compatible assembly/disassambly
- > Exchangeable insulation sleeves available (PEEK/PEI)
- > Customizable to meet your requirements
- > Helium-leakage rate < 10^{-7} mbar l s⁻¹ (at 250 °C)





CompreDriveControl – Automate your measurements



Many devices – connected by one program

- Interconnection of thermocontroller, analyzer software and CompreDrive
- > Definition of individual, automated procedures

Comprehensive control and software interface

- > Full manual control over force and temperature setpoints
- > Extensive data recording and export options
- > Ability to implement and record additonal devices via plugins
- > Set force or pressure using built-in cell definitions





Automated measurement routines

- Procedures defined by setpoint lists
- > Customizable end-conditons for setpoint steps
- > Synchronization with analyzers by trigger signals

HTTP-based software interface

- Control the CompreDrive with HTTP request from external software
- Remote setpoints directly from analyzer software (if supported)

System requirements (recommended)

rhd III instruments

flexible cell solutions

| Operating system: | Windows 10 |
|-------------------|---------------------------|
| RAM: | 2 GB |
| Hard disk space: | 200 MB |
| Ports: | USB 2.0 (for CompreDrive) |

Compre Dr ve


Specifications:

| Maximum force: | 75 kN |
|-------------------------------|--|
| Force precision: | ±0.1 kN |
| Maximum pressure: | 1.7 GPa (Ø 6 mm)/663 MPa (Ø 12 mm) |
| Temperature compliance range: | -40 °C \leftrightarrow +250 °C** **depending on sleeve material |
| Temperature accuracy: | ±0.1 °C |
| Materials in sample contact: | Hard-metal, AlOx-ceramic/PEEK/PEI |
| Sample diameter: | 6 or 12 mm |
| Sample height: | ≤ 10 mm |
| Drive type: | Servo-driven lifting spindle |
| Power supply: | 1-phase AC 110-230 V, 50/60 Hz |
| Maximum total cell height: | 149 mm |
| Dimensions: | 167 x 81 x 58 cm |
| Total weight: | 180 kg |



rhd instruments

flexible cell solutions





For questions, quotes and orders please contact us:



rhd instruments GmbH & Co. KG Otto-Hesse-Strasse 19 / T3 64293 Darmstadt, Germany

info@rhd-instruments.de www.rhd-instruments.de +49 6151 870 7187



08/2020